

ACCESSORIES

Batteries & Transformers for Bionet Veterinary Monitors


www.BionetUS.com
714-734-1760
support@BionetUs.com

COMPATIBLE DEVICES




BM Rechargeable battery: LEAD-ACID
Item Code: BM-BAT-1

BM3: D1C0400001 ~ D1E08xxxxx(04.2002 ~ 08.2004)
BM3Vet: D1C0400001 ~ D1E08xxxxx(04.2002 ~ 08.2004)
BM3VetPlus: D2D0200001 ~ D2D07xxxxx(2003.02 ~ 07.2003), BMP0308 ~ BMP0408(BMP0308 ~ BMP0408)
BM3Plus: D2D0200001 ~ D2D07xxxxx(2003.02 ~ 07.2003), BMP0308 ~ BMP0408(BMP0308 ~ BMP0408)



Rechargeable battery: Ni-MH (2line, each line is 5cell, 12V, 1800mA/h)
Item Code: BM-BAT-2

BM3: D1E0900001 ~ D1F11xxxxx(09.2004 ~ 11.2005)
BM3Vet: D1E0900001 ~ D1F11xxxxx(09.2004 ~ 11.2005)
BM3VetPlus: BMP0409 ~ BMP0509, D2F1000001 ~ D2F11xxxxx(09.2004 ~ 11.2005)
BM3Plus: BMP0409 ~ BMP0509, D2F1000001 ~ D2F11xxxxx(09.2004 ~ 11.2005)



Rechargeable battery: Ni-MH (10 cell, 1line, 12V, 2200mA/h)
Item Code: BM-BAT-3

BM3: D1F1200001 ~ D1G03xxxxx(12.2005 ~ 03.2006), D1G04 ~ D1H0300001 - 57(04.2006 ~ 03.2007), D1H0300058 - 110 ~ D1I0300150(03.2007 ~ 03.2008)
BM3VetPlus: D2F1200001 ~ D2H02xxxxx(12.2005 ~ 02.2007)
BM3Vet: D1F12 ~ D1I03xxxxx(12.2005 ~ 03.2008)
BM3Plus: D2F1200001 ~ D2H02xxxxx(12.2005 ~ 02.2007)



BM Rechargeable battery: Li-ION, 3cell, 1line, 10.8V, 2200mA/h
Item Code: BM-BAT-4

BM3: D1I0400001 ~ UP TO DATE (04.2008 ~ UP TO DATE)
BM3Vet: D1I0400001 ~ UP TO DATE (04.2008 ~ UP TO DATE)
BM3: DDL0700001 ~ UP TO DATE (07.2011 ~ UP TO DATE)
BM3Vet Next: DDL0700001 ~ UP TO DATE (07.2011 ~ UP TO DATE)
BM3VetPlus: D2H0300001 ~ D2I07XXXXXX(03.2007 ~ 07.2008), D2I1200001 ~ D2M12xxxxx(12.2008 ~ 12.2012)
BM3Plus: D2H0300001 ~ D2I07xxxxx(03.2007 ~ 07.2008), D2I1200001 ~ D2M12xxxxx(12.2008 ~ 2012)
BM5Vet: D4G0900001 ~ D4I10xxxxx(09.2006 ~ 10.2008)
BM5Vet: D7J0400001 ~ D7K10xxxxx(04.2009 ~ 10.2010)
BM5Vet: D8K1100001 ~ UP TO DATE(11.2010 ~ UP TO DATE)
BM3Vet Touch: All units



Rechargeable battery: Li-ION (EtCo2) 6 Cell, 2 Line, 11.1V, 4400mA/h
Item Code: BM-BAT-5

BM5Vet: D4G0900001 ~ D4I10xxxxx(09.2006 ~ 10.2008)



Rechargeable battery: Li-Ion, 6Cell, 1Line, 10.8V, 4300mAh/46.44Wh
Item Code: BM-BAT-6

BM7Vet: DEP0400096 ~ UP TO DATE
BM7Vet Pro



Rechargeable Li-ion battery, 3.6V, 4300mAh/15.48Wh
Item Code: OXW-BAT

Oxy9Wave Vet: O2P040090 ~ UP TO DATE



BM Rechargeable battery: Li-ION 3 Cell, 1 Line, 11.1V, 2200mA/h
Item Code: BM1-BAT

BM1: DBJ1000001 ~ UP TO DATE(10.2009 ~ UP TO DATE)
BM1 (Suntech): DCL0800001 ~ UP TO DATE(08.2011 ~ UP TO DATE)



Optional Rechargeable Battery: Li-ion, 6Cell 2Lines, 10.8V, 4800mAh
Item Code: BM-BAT-7

BM7Vet: DEP0400096 ~ UP TO DATE
BM7Vet Pro

COMPATIBLE DEVICES



INPUT: 100~240V, 50-60 Hz, 1.5A
OUTPUT: +18 VDC, 2.8A
(Doesn't include power cable (PC-A(V)))
Item Code: TRANS

BM3, BM3PLUS, BM3VET, BM3VET TOUCH, BM5VET, BM7VET, BM PRO SERIES (BM3VET PRO / BM5VET PRO / BM7VET PRO)



INPUT: 100~240V, 50-60 Hz, 1.0A
OUTPUT: +15 VDC, 2.0A
(Doesn't include power cable (PC-A(V)))
Item Code: BM1 TRANS

BM1VET



INPUT: 100~240V, 50-60 Hz, 0.4A
OUTPUT: +5 VDC, 2.0A
(Doesn't include power cable (PC-A-OX))
Item Code: OXW-TRANS

OXY9WAVE VET



110V power cable
Item Code: PC-A(V)

TRANS, BM1VET TRANS



110V power cable
Item Code: PC-A-OXW

OXW-TRANS



INPUT: 110~240V, 47-63 Hz, 0.48A
OUTPUT: +9 VDC, 1.5A
(Doesn't include power cable (PC-A-OX))
Item Code: OXP-TRANS

Oxy9Vet Plus



BM Rechargeable battery: Li-ION 3 Line, 10.8V, 2150mA/h
Item Code: BM-BAT-3&5

BM3Vet Pro, BM5Vet Pro



BM Rechargeable battery: Li-ION 1 Line, 3 Cell, 10.8V, 3250mA/h
Item Code: BM-BAT-3&5(N)

BM3Vet Elite, BM5Vet Elite



BM Rechargeable battery: Li-ION 1 Line, 6 Cell, 10.8V, 6500mA/h
Item Code: BM-BAT-7(N)

BM7Vet Elite