

Configurations to fit hospital of all sizes and operations

Model	Super Air C1 Series			Super Air C2 Series			Super Air C4
	 Air-C1-S (Single)	 Air-C1-D (Single w/stand)	 Air-C1-T (Twin)	 Air-C2-S (Single)	 Air-C2-D (Single w/stand)	 Air-C2-T (Twin)	 Air-C4-T
Chamber	1 Chamber		2 Chambers	1 Chamber	2 Chambers		3 Chambers
Chamber Internal Dimension	27.2"W x 19.5"D x 20.5"H		Upper & lower, each 27.2"W x 19.5"D x 20.5"H	42.0"W x 24.6"D x 24.0"H (Dividable)	Upper & lower, each 42.0"W x 24.6"D x 24.0"H (Dividable)		Upper: 2 Chambers 21.5"W x 28.0"D x 20.4"H Lower: 1 Chamber 47.9"W x 28.0"D x 28.3"H (Dividable)
Controller	1 controller		2 controllers	1 controller	2 controllers		4 controllers
Maximum patient	1 patient		2 patients	1~2 patients	2~4 patients		3~4 patients
Stand (Stroage)	No	Yes	No	No	Yes	No	No
External Dimension	34.8"W 30.3"D 35.8"H	34.8"W 30.3"D 67.6"H (including Stand)	34.8"W 30.3"D 67.6"H	49.3"W 34.8"D 39.0"H	49.3"W 34.8"D 74.2"H (including Stand)	49.3"W 34.8"D 74.2"H	55.2"W 34.8"D 75.4"H